



An ISO 13485:2016 Certified company

Govt. of India CDSCO License : MFG/MD/2023/000918

PULSELITH PNEUMATIC LITHOTRIPTER

PART NO. AHCR-PLT-001-001

- Inlet supply : 6Kgs/ cms Square (max)
- Inlet supply : 4Kgs/ cms Square(min)
- Input Gas : Oxygen/Compressed Air
- Output Pressure : 2.5to4.5Kgs cms Square
- Mains Supply Voltage : 180 Vac - 230 Vac 50 Hz
- Power Consumption : 40 Watts
- Frequency : Single/Multiple/Power
- Projectile Impact Frequency : 1 Hz - 12 Hz and 14 Hz
- Probes-0.8mm, 1.0mm,1.2mm, 1.5mm, 2.0mm, 2.5mm



ACCESSORIES OF THE PNEUMATIC LITHOTRIPTER



CARRY BAG WITH COMPLETE SYSTEM



AHCR-PLT-001-001

Office & Factory : B-121, 3rd Floor, DDA Shed, Okhla Industrial Area,
Phase-1, New Delhi-110 020 INDIA,  www.pulselith.com



An ISO 13485:2016 Certified company

Govt. of India CDSCO License : MFG/MD/2023/000918

PULSELITH SILENT AIR COMPRESSOR

AHCR-COM-001-001

- Oil Based Sealed Compressor
- ½ HP Compressor
- Duty Cycle 50 HZ
- Tank capacity 5 L
- Maximum flow rate 120 L per minute
- Pressure adjustable from 1 to 7 bars
- Wattage 350 watts
- The entire unit is on castors



AHCR-COM-001-001

Office & Factory : B-121, 3rd Floor, DDA Shed, Okhla Industrial Area,
Phase-1, New Delhi-110 020 INDIA,  www.pulselith.com