



ADVANCED
HEALTH CARE RESOURCES
Quality Equipment • Premium Services • Competitive Pricing

An ISO 13485:2016 Certified company

Govt. of India CDSCO License : MFG/MD/2023/000918

LAPAROSCOPIC CO2 INSUFFLATOR 30 L

AHCR-CIS-001-001



- 30L Insufflators Precise control of pressure/flow Parameters.
- Greater accuracy & measurement of Co2
- Insufflation Pressure :- 3-25 mm Hg
- Gas Connection :- UNF 7/16"
- Input Pressure :- 3 Bar min.
- Output Connection :- 9 - 10 mm ID Tube
- Flow Rate Setting :- 1 to 30 L/ Min



AHCR-CIS-001-001

Office & Factory : B-121, 3rd Floor, DDA Shed, Okhla Industrial Area,
Phase-1, New Delhi-110 020 INDIA,  www.pulselith.com