

CYSTOSCOPY SET

ADVANCED HEALTH CARE RESOURCES



ADVANCED
HEALTH CARE RESOURCES
Quality Equipment • Premium Services • Competitive Pricing

Scope Ø 4.0mm, 30° viewing angle, autoclavable




Single Channel Bridge




Double Channel Bridge



Cystoscopy sheath 17 Fr, 19 Fr, 20, 21 Fr, 23.5 Fr & 25 Fr. with 2 fixed stopcock, with obturator



*Biopsy forceps (5Fr & 7Fr)
Biopsy forceps Flexible and Semi Rigid*



*Cystoscopy Alligator Jaw (7Fr)
Flexible and Semi Rigid*